

Lecture #2: Enzyme Purification and Assay

(1) Initial velocity measurements

- Effect of enzyme concentration
- Effect of temperature

(2) Assay types

- Discontinuous
- Continuous
- Quantitative assays
- Coupled reactions

(3) Enzyme activity measurements

- Units
- Katal

(4) Turnover number and properties

(5) Purification and assessment

- Purification table


(6) Methods for measurements

- Continuous systems
- Immobilized enzymes
- Stopped flow
- Quench flow
- Multi-plate readers



Lecture #2: Enzyme Purification and Assay

How to Perform Kinetic Measurements



1. Dead time

- time after mixing the sample and before the measurement begins

2. Initial slope

- Linear portion of the time-course plot for data collection

3. Post collection

- Time after initial velocity during which the reaction reaches equilibrium

1. Initial velocity measurements




Figure 6-10
Lehninger Principles of Biochemistry, Seventh Edition
© 2017 W. H. Freeman and Company

• a typical enzyme reaction is reduced as the reaction proceeds because:



- denaturation of enzyme
- product inhibition
- decrease of enzyme saturation
- inactivation of coenzyme
- increase of reverse reaction

• critical that **velocity** be determined at the very beginning of the reaction



• initial velocity is the slope of the curve depicting S or P changes with time

$$v_o = d[P]/dt = -d[S]/dt$$

• in practice: $v_o = \Delta[P]/\Delta t = -\Delta[S]/\Delta t$ **provided that** these parameters are measured at the **beginning** of the **reaction** where the rate of S consumption and P formation is linear with time



- By choosing an arbitrary incubation time, the incubation time itself now becomes a variable
- Can be avoided by measuring initial velocities




Core Topics in
Biochem.,
Stenesh (1993)

2. Assay types

i. Discontinuous assay - fixed time point sampling of kinetic data


- Radiometric: measure incorporation of radioactivity into substrates
- Chromatographic: measure product formation by separating the reaction mixture into components by chromatography (HPLC, FPLC, TLC)

•the rate of [P] formation must be linear over the incubation time chosen



ii. Continuous assay - reaction is analyzed continuously using a monitoring technique

- Spectrophotometric, calorimetric, chemiluminescent
- for the continuous assay one must use the initial part of the curve for dP/dt
- one must use [E] that falls within the range in which P formation varies linearly with [E] for all substrate concentrations



Two types of quantitative enzyme assays

(1) Measure the amount of E in the sample

- E is present at **much lower** concentration than S
- E is present in **limiting amounts** and S is present **in excess**
- assay conditions correspond to levelled off portions of v_o versus $[S]$ curves

thus velocity $\propto [E]$





Figure 6-11
Lehninger Principles of Biochemistry, Seventh Edition
© 2017 W. H. Freeman and Company




Zero order reaction

$$V_o = [E] * k * [S]^0$$

(2) Measure the amount of S present ([S])

- [E] is high enough to convert all $S \rightarrow P$
- S is limiting reagent
- assay conditions correspond to those of initial part of v_0 versus [S] curves



First order reaction





Figure 6-11
Lehninger Principles of Biochemistry, Seventh Edition
© 2017 W. H. Freeman and Company

$$V_0 = [E] * k * [S]^1$$

(3) Coupled Reactions

- some reactions difficult to assay because of **absence of a readily measurable** compound or property
- couple one reaction to another which produces a measurable product



- if conditions for E_2 are compatible with those of E_1 then both stages can be measured **simultaneously**
- if not, then reaction with E_1 is run first and stopped (boiling, change pH, etc.)
- second reaction (with E_2 in excess) is added and reaction proceeds until all of B is converted to C
- for such "coupled reactions" E_2 must be in **sufficient concentration** so that all of B is converted to C
 - Rate-limiting step in the linked assay is the **action of the first enzyme**

- the relative amounts of E_1 and A depend on the nature of assay
 - if E_1 is being measured then A should be in **high concentration**
 - if A is being measured then A should be in **low concentration**

3. Enzyme units of activity

- enzyme solutions described in terms of arbitrary units because enzymes are **active at very low concentrations** and purity is frequently unknown



Unit (definition)

- one international unit = amount of enzyme catalyzing the transformation of **1 μmol of substrate per min** under defined conditions
 - based on measurements of **initial velocity**

-greater the **purity** of an enzyme, the greater the fraction of total protein that is ***actually enzyme protein***

-relationship defined as "**specific activity**"

-**specific activity** - number of enzyme units per amount (mg) of total protein= **U/mg protein**

-**total activity** = U/mg protein x total mg protein
= U/mL of fraction x total vol (mL) of fraction



I.U.B. proposed **katal** (1972) \equiv **SI Unit of activity**

katal is amount of enzymatic activity that catalyzes the transformation of **1 mole of substrate per second** under defined conditions

$$1 \text{ U} = 10^{-6} \text{ moles/60s} = 16.7 \times 10^{-9} \text{ mol/s}$$

$$\text{- } 1 \text{ U} = 16.7 \text{ nkatal}$$

$$\text{- } 1 \text{ nkatal} = 0.06 \text{ U}$$

$$\text{- } 1 \text{ katal} = 6 \times 10^7 \text{ U}$$

-concentration of enzyme in terms of katal is "**molar activity**"

-defined as the number of *katal*s per mole of enzyme

4. Turnover number (TN, k_{cat})

-number of moles of substrate transformed into product per unit time per mole of enzyme under optimal conditions

-**molecular activity** = TN (per mole of enzyme)

-**catalytic centre activity** = TN (per mole of active sites)

$$\begin{aligned}\bullet \text{TN} &= V_{\text{max}}/[E_T] = \mu\text{mol}(\text{S} \rightarrow \text{P}) \times \text{min}^{-1} \times \text{mL}^{-1} / (\mu\text{mol enzyme or cat site} \times \text{mL}^{-1}) \\ &= \text{min}^{-1} \text{ or } \text{t}^{-1}\end{aligned}$$

-typical values fall in the range of 50 - 10^7 min^{-1}

TABLE 6-7 Turnover Number, k_{cat} , of Some Enzymes

Enzyme	Substrate	k_{cat} (s^{-1})
Catalase	H_2O_2	40,000,000
Carbonic anhydrase	HCO_3^-	400,000
Acetylcholinesterase	Acetylcholine	14,000
β -Lactamase	Benzylpenicillin	2,000
Fumarase	Fumarate	800
RecA protein (an ATPase)	ATP	0.5

5. Purification and assessment

-purification is the process by which an enzyme is separated from other cell components




Three categories of purification methods:

i. **Precipitation** \Rightarrow differential solubility

- salting out with $(\text{NH}_4)_2\text{SO}_4$ most common.
- salts neutralize exterior charge of protein allowing aggregation.
- isoelectric precipitation uses same principle

ii. **Chromatography**

- various types, including
 - gel filtration
 - ion-exchange
 - hydrophobic
 - adsorption



- affinity: taken over with recombinant enzymes today!
 - His tags, Flag tags, GST-tags

iii. **Electrophoresis**


- native or denaturing

-presentation of purification data is vital for analysis of purity and recovery

-aim is to calculate **yield** and **purification factor**

- **need 3 measurements:**

- i. volume of sample (mL)
- ii. protein content (mg/mL)
- iii. enzyme activity (units/mL)




- **all other quantities derived by calculation:**

$$\text{sample vol (mL)} \times \text{protein concn (mg/mL)} = \text{total protein (mg)}$$

$$\text{sample vol (mL)} \times \text{activity (units/mL)} = \text{total activity (units)}$$

$$\frac{\text{total activity (units)}}{\text{total protein (mg)}} = \text{specific activity (units/mg)}$$

Figure 3-18b
Lehninger Principles of Biochemistry, Sixth Edition
© 2013 W. H. Freeman and Company



I found that the air kills it.
Room temperature, too.

Now it dies in the cold room.
I've got the enzyme purification blues.

- the higher the yield, the greater the recovery, the fewer times you'll need to get in the cold room and purify the enzyme
- the higher the purification factor, the purer the enzyme
 - the purification factor still doesn't give information about **how pure** it really is!

TABLE 3-5 A Purification Table for a Hypothetical Enzyme

Procedure or step	Fraction volume (mL)	Total protein (mg)	Activity (units)	Specific activity (units/mg)
1. Crude cellular extract	1,400	10,000	100,000	10
2. Precipitation with ammonium sulfate	280	3,000	96,000	32
3. Ion-exchange chromatography	90	400	80,000	200
4. Size-exclusion chromatography	80	100	60,000	600
5. Affinity chromatography	6	3	45,000	15,000

Note: All data represent the status of the sample *after* the designated procedure has been carried out. Activity and specific activity are defined on page 95.

Table 3-5

Lehninger Principles of Biochemistry, Sixth Edition

© 2013 W. H. Freeman and Company

What properties of any enzyme should be determined?

(i) molecular weight

- electrospray mass spectrometry
- gel filtration
- analytical ultracentrifugation

(ii) subunit structure

- SDS-PAGE (\pm DTT)

(iii)



- amino acid composition
- 6M HCl, 105°C, 24 h

(iv)

- amino acid sequence
- Edman degradation
- N-terminal sequencing
- MS/MS mass spectrometry

(v) -three dimensional structure


- X-ray crystallography
- NMR



6. Methods for Enzyme Activity Measurements

(a) Experimental approaches

- i. **spectrophotometric assays**- accurate and convenient
 - fluorescence or absorbance changes
- ii. **Release or uptake of H^+ (ΔpH)**- directly in buffered or weakly buffered solutions with a glass electrode
 - restricted pH range (enzyme stability)
 - indicator that changes its absorbance with protonation state
 - pH stat technique**-titrate the reaction mixture with either acid or base to keep the pH constant whilst recording the rate of addition
 - ion-selective electrodes or gas electrodes




iii. **Automated assay procedures**

- automation of E assay can allow **large numbers of samples** to be processed rapidly and efficiently
- many involve the determination of P formation after a fixed time

- (1) **continuous flow systems**-uses multi-channel pumps to mix reactants and determine P formation after a fixed time

On-line enzyme MS assay – System 1






(2) **immobilized enzymes**-enzymes are attached to a matrix (bead or surface)


eg., trypsin-bead columns

(3) **stopped-flow apparatus**

-enzyme mixed with substrate using a flow cell and detected with spectrometer


-allows rapid kinetic measurements (msec time scale)





- Stopped-flow kinetic trace for the binding of a substrate to the enzyme
- Substrate binding causes quenching of the enzyme Trp fluorescence

4) **quenched-flow apparatus**-enzyme is mixed with substrate and quenched with acid or base in a second mixer



5) **multi-plate readers**-wavelength selection and temperature control coupled with rapid sample throughput

McMaster HTS Lab



<http://www.cmcbmcmaster.ca/high-throughput-screening-lab/>

